

KRT18

Mouse Anti-Human Keratin 18 (Clone C-04) mAb

Catalog No.	MON3034	Quantity:	100 ug
Alternate Names:	Keratin, type I cytoskeletal 18, Cell proliferation-inducing gene 46 protein, Cytokeratin-18, CK-18, Keratin-18, K18		
Description:	MON 3034 reacts with Cytokeratin 18 (45 kDa), a member of intermediate filaments subfamily represented in epithelial tissues.		
Concentration:	1 mg/ml		
Uniprot ID:	P05783		
Specificity:	Keratin 18		
Host:	Mouse		
Immunogen:	Cytoskeleton preparation of epidermal carcinoma cell line A431		
Isotype:	IgG1		
Clone:	C-04		
Formulation:	Purified antibody in PBS containing 15 mM sodium azide		
Applications:	WB, IP, ICC, IHC-P, FC		
Storage & Stability:	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.