

Mouse Anti-Human Basic Keratin Type II (Clone AE-3) mAb

Catalog No.	MON3038	Quantity:	250 ug
Description:	MON 3037 recognizes all known basic (Type II) keratins. It's broadly reactive and stains positively almost all epithelia and their neoplasms. This antibody together with MON 3037 has been used for studying the immuno-localization of epidermal keratins, tissue distribution and epithelia development. AE-1 and AE-3 have also been used extensively, individually or mixed, for distinguishing carcinoma from noncarcinoma		
Concentration:	1 mg/ml		
Specificity:	Epithelial Keratin Type II		
Host:	Mouse		
Immunogen:	Human epidermal cells		
Isotype:	IgG1		
Clone:	AE-3		
Formulation:	Purified antibody in PBS containing 0.1% sodium azide		
Cross-Reactivity:	Reacts with keratin from rabbit, mouse, rat, bovine and chicken.		
Applications:	Western blot: suggested dilution 1:100 IHC: suggested dilution 1:150 - 1:200		
Storage & Stability:	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

