

## H1FO

### Mouse Anti-Human Lewis A (Clone 7LE) mAb

<b>Catalog No.</b>	MON3051	<b>Quantity:</b>	100 µg
<b>Description:</b>	MON3051 reacts with Lewis A blood group antigen, a carbohydrate determinant present on both glycolipids and glycoproteins, expressed a.o. in colonic and kidney epithelial cells. Lewis A may be useful for detection of gastrointestinal, pancreatic and colorectal carcinomas.		
<b>Concentration:</b>	0.1 mg/ml		
<b>Specificity:</b>	Lewis A		
<b>Host:</b>	Mouse		
<b>Immunogen:</b>	Ovarian cyst mucin		
<b>Isotype:</b>	IgG1 kappa		
<b>Clone:</b>	7LE		
<b>Formulation:</b>	Purified antibody in PBS with 0.02% sodium azide.		
<b>Applications:</b>	IHC, ELISA, FC, WB, IF		
<b>Storage &amp; Stability:</b>	Store antibody at 2-8°C until expiration date.		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

