

Mouse Anti-Human Cancer Antigen 19-9 mAb (Clone 121SLE)

Catalog No.	MON3060	Quantity:	100 µg
Alternate Names:	CA19-9, Sialyl Lewis A		
Description:	MON 3060 recognizes Sialyl Lewis A containing glycolipids. High levels of CA19-9, indicate malignancies of the pancreas, colorectum, lung, liver, and ovary.		
Concentration:	1 mg/ml		
Specificity:	CA19.9		
Host:	Mouse		
Immunogen:	Immunoprecipitate obtained after immunodiffusion of MAb 19-9 and mucins isolated from an ovarian cyst of a OLe(a+b-) patient		
Isotype:	IgM		
Clone:	121SLE		
Formulation:	Purified antibody in Tris buffered saline (TBS) with 15 mM NaN ₃		
Applications:	ELISA, FC, IHC-P		
Storage & Stability:	Store undiluted antibody at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.