

## Mouse Anti-Human Keratin 6, (Clone Ks6.KA12) mAb

<b>Catalog No.</b>	MON3070	<b>Quantity:</b>	1.0 mL
<b>Alternate Names:</b>	Keratin, type II cytoskeletal 6, KRT6		
<b>Description:</b>	MON3070 represents an excellent marker for non-keratinized squamous epithelia and proliferating cells of epidermis (e.g. within psoriatic lesions). Reacts with 56 kDa polypeptide (cytokeratin 6) of human squamous epithelia.		
<b>Specificity:</b>	Keratin 6		
<b>Host:</b>	Mouse		
<b>Immunogen:</b>	Keratin K6 of human callus cytoskeletal preparation		
<b>Isotype:</b>	IgG1		
<b>Clone:</b>	Ks6.KA12		
<b>Formulation:</b>	PBS containing 0.5% BSA and 0.09% sodium azide		
<b>Applications:</b>	WB, IHC		
<b>Application Notes:</b>	This antibody is suitable for discrimination of keratinizing and non-keratinizing squamous cell carcinoma versus e.g. poorly differentiated adenocarcinoma.		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.