

Mouse Anti-Human Myelin (Clone 2B5) mAb

Catalog No.	MON4004	Quantity:	1.0 mL
Description:	MON 4004 recognizes in paraffin sections only myelin in the peripheral nervous system (PNS). In frozen sections it also reacts with myelin in the central nervous system (CNS). Is useful for studying myelination and demyelination processes in the PNS, and remyelination processes in the CNS as can be observed after trauma in multiple sclerosis. Reacts with a 43 kD protein in immunoblotting		
Specificity:	Myelin		
Immunogen:	Human myelin isolated from peripheral nervous system		
Isotype:	Mouse IgM		
Clone:	2B5		
Formulation:	Ascites diluted in PBS containing 0.1% BSA and 0.09% sodium azide		
Applications:	IHC: 1:10 WB		
Storage & Stability:	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

