

Mouse Anti-Human Myosin (Clone 414) mAb

Catalog No.	MON4005	Quantity:	1.0 mL
Description:	MON4005 reacts with myosin-containing elements in most cell types. Antigen localization: cytoplasm In Western blot experiments the antibody reacts with a 75kD part of the heavy chain (epitope location on a myosin-common site-chain), no reaction with light chains. It is not reactive with alpha-cardiac myosin.		
Specificity:	Myosin		
Immunogen:	Chicken skeletal muscle preparation		
Isotype:	Mouse IgM		
Clone:	414		
Formulation:	Culture supernatant containing 0.09% sodium azide		
Applications:	IHC: suggested dilution 1:5 - 1:10 WB, IF, ICC		
Storage & Stability:	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.