

Mouse Anti-Human Beta-Tubulin (Clone TU-06) mAb

Catalog No.	MON4010	Quantity:	100 ug
Description:	MON4010 reacts with β -tubulin of various species including that of plants. It reacts with the N-terminal structural domain of β -tubulin. This antibody reacts with all isoforms of β -tubulin. Clone TU-06 reacts with all charge variants of tubulin.		
Concentration:	1.0 mg/ml		
Specificity:	Beta-Tubulin from various species including plants		
Immunogen:	Beta-subunits of porcine brain tubulin		
Isotype:	Mouse IgM		
Clone:	TU-06		
Formulation:	PBS containing 15 mM sodium azide.		
Cross-Reactivity:	Reacts with beta-tubulin in human, mouse, rat, hamster, pig, cow, chicken, carp as well as in plants, yeasts and protozoa		
Applications:	IHC-P, ICC, WB, IP		
Storage & Stability:	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.