

## Mouse Anti-Human Alpha-Tubulin (Clone TU-16) mAb

<b>Catalog No.</b>	MON4012	<b>Quantity:</b>	100 ug
<b>Description:</b>	MON 4012 binds tubulin under both denaturing and non-denaturing conditions. The antibody reacts with tubulins of various species.		
<b>Concentration:</b>	1.0 mg/ml		
<b>Specificity:</b>	Alpha-Tubulin		
<b>Immunogen:</b>	Porcine brain microtubule protein MTP-1		
<b>Isotype:</b>	Mouse IgM		
<b>Clone:</b>	TU-16		
<b>Formulation:</b>	PBS containing 15 mM sodium azide		
<b>Cross-Reactivity:</b>	Cross-reactivity was observed in human, mouse, rat, hamster, pig, cow, chicken as well as in plants and protozoa.		
<b>Applications:</b>	ICC, IHC-P, IP, WB		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 year.		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**