

TUBA3

Mouse Anti-Human Tubulin-alpha (Clone TU-02) mAb

Catalog No.	MON4018	Quantity:	100 µg
Description:	MON4018 reacts with the N-terminal domain of alpha-tubulin. Tubulin isotypes have been identified with tissue specific expression. Immunocytochemical studies using several mAbs revealed remarkable heterogeneity of tubulin within various nervous tissues. It reacts with tubulin in fixed tissues only, not in unfixed or live tissues or cells. Interphase microtubules are also stained by TU-02 in fixed tissues.		
Concentration:	0.1 mg/ml		
Specificity:	N-terminal domain of Alpha-Tubulin		
Immunogen:	Microtubular protein from pig brain		
Isotype:	Mouse IgM		
Clone:	TU-02		
Formulation:	Purified antibody in PBS with 0.02% sodium azide.		
Applications:	ELISA: detecting antibody at 0.001 - 100 µg/ml IHC-P: 1-2 µg/ml, 30 min at RT		
Storage & Stability:	Store antibody at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

