

## LAT

### Mouse Anti-Human LAT (Clone LAT-01) mAb

<b>Catalog No.</b>	MON4023	<b>Quantity:</b>	100 µg
<b>Alternate Names:</b>	Linker for activation of T-cells family member 1, 36 kDa phospho-tyrosine adapter protein, pp36, p36-38		
<b>Description:</b>	MON 4023 reacts with an intracellular epitope of LAT, a 36-38 kDa transmembrane adaptor protein expressed by T cells, pre-B cells, NK cells, mast cells and platelets.		
<b>Concentration:</b>	1 mg/ml		
<b>Gene ID:</b>	O43561		
<b>Specificity:</b>	LAT		
<b>Immunogen:</b>	Bacterially produced recombinant polypeptide corresponding to the entire cytoplasmic domain of human LAT		
<b>Isotype:</b>	Mouse IgG1		
<b>Clone:</b>	LAT-01		
<b>Formulation:</b>	PBS containing 15 mM sodium azide, pH 7.4		
<b>Applications:</b>	WB, FC, IP		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 year.		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**