

## AFP

### Mouse Anti-Human Alpha-fetoprotein (Clone C3) mAb

<b>Catalog No.</b>	MON4036	<b>Quantity:</b>	100 µg
<b>Alternate Names:</b>	Alpha-fetoprotein, Alpha-1-fetoprotein, Alpha-fetoglobulin		
<b>Description:</b>	Alpha-fetoprotein is one of the major serum proteins in the early life of mammals. AFP is an oncodevelopmental gene product which is expressed at high levels in the embryonic yolk sac and fetal liver. AFP production is reactivated in the adult during liver regeneration and hepatocarcinogenesis.		
<b>Concentration:</b>	0.1 mg/ml		
<b>UniProt ID:</b>	P02771		
<b>Specificity:</b>	Alpha Fetoprotein (AFP)		
<b>Host:</b>	Mouse		
<b>Immunogen:</b>	Alpha fetoprotein isolated from serum of hepatoma patients		
<b>Isotype:</b>	IgG2a		
<b>Clone:</b>	C3		
<b>Formulation:</b>	Sterile-filtered purified antibody in PBS with 0.1% BSA and 0.02% sodium azide		
<b>Applications:</b>	IHC, ELISA, WB		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.