

AFP

Mouse Anti-Human Alpha-fetoprotein (Clone D10) mAb

Catalog No.	MON4037	Quantity:	100 µg
Alternate Names:	Alpha-fetoprotein, Alpha-1-fetoprotein, Alpha-fetoglobulin		
Description:	Alpha-fetoprotein is one of the major serum proteins in the early life of mammals. AFP is an oncodevelopmental gene product which is expressed at high levels in the embryonic yolk sac and fetal liver. AFP production is reactivated in the adult during liver regeneration and hepatocarcinogenesis.		
Concentration:	0.1 mg/ml		
UniProt ID:	P02771		
Specificity:	Alpha Fetoprotein (AFP)		
Host:	Mouse		
Immunogen:	Alpha fetoprotein isolated from serum of hepatoma patients		
Isotype:	IgG1		
Clone:	D10		
Formulation:	Sterile-filtered purified antibody in PBS with 0.1% BSA and 0.02% sodium azide		
Applications:	IHC, ELISA, WB		
Storage & Stability:	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.