

ITGAE

Mouse Anti-Human CD103 (Clone LF61) mAb

Catalog No.	MON4055	Quantity:	100 µg
Alternate Names:	Integrin alpha-E, HML-1 antigen, Integrin alpha-IEL, Mucosal lymphocyte 1 antigen		
Description:	Recognizes human CD103, the cell surface Integrin alpha-E antigen, a glycoprotein expressed by approximately 1% of peripheral blood lymphocytes, activated T lymphocytes and by hairy cell leukaemia cells. The antigen is also expressed by intraepithelial lymphocytes.		
Concentration:	1 mg/ml		
UniProt ID (target):	P38570		
Specificity:	CD103		
Host:	Mouse		
Immunogen:	Hairy cell leukemia cells		
Isotype:	IgG ₁		
Clone:	LF61		
Formulation:	PBS containing 0.09% sodium azide.		
Applications:	FC, IHC-F, IP		
Storage & Stability:	Store at 2-8°C for 1 month or in undiluted aliquots at -20°C to -80°C up to 1 year. Avoid repeated freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

