

Mouse Anti-Human Lambda Light Chain (Clone 48) mAb

Catalog No.	MON5002-1	Quantity:	1 ml
	MON5002-5		5 ml

Description: This antibody will react with free lambda chains and also with lambda chains attached to heavy chains. In Ouchterlony immunodiffusion, no reaction with Bence Jones lambda or Ig lambda. No cross reactivity with kappa chains.

Concentration: 20 ug/ml

Specificity: Human lambda light chain

Host: Mouse

Immunogen: Lambda Bence Jones protein

Isotype: IgG₁

Clone: Clone 48

Formulation: PBS containing 10 mM sodium azide and 1% FCS.

Applications: IHC

Storage & Stability: Store at 2-8°C for 1 month or in undiluted aliquots at -20°C to -80°C up to 1 year. **Avoid repeated freeze/thaw cycles.**

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

