

Mouse Anti-Beta-Galactosidase, *E. coli* (Clone BG-02) mAb

Catalog No.	MON5013	Quantity:	100 ug
Description:	This antibody recognizes a non-active center epitope of beta-galactosidase from <i>E. Coli</i> . Cross reactivity with LAC-Zalpha fragment is unknown. The antibody reacts with lambda gt11 fusion proteins.		
Concentration:	1 mg/ml		
Specificity:	Beta-Galactosidase from <i>E. coli</i>		
Host:	Mouse		
Immunogen:	Beta-Galactosidase purified from <i>E. coli</i>		
Isotype:	Murine IgG1		
Clone:	BG-02		
Formulation:	Ascites diluted in PBS containing 0.1% sodium azide		
Applications:	WB, ICC, IF, IHC		
Application Notes:	The antibody can be used in molecular biological tests, immuno-enzymological studies and tests: galactosidase-antigalactosidase (GAG) technique. Also applicable for isolation of lambda gt11 fusion proteins.		
Storage & Stability:	Store at 2-8°C for 1 month or in undiluted aliquots at -20°C to -80°C up to 1 year. Avoid repeated freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

