

ALB

Mouse Anti-Human Serum Albumin (Clone AL-01) mAb

Catalog No.	MON5024	Quantity:	100 ug
Alternate Names:	HSA		
Description:	This antibody is highly specific for human albumin. No reactivity with other serum proteins. The antibody does not cross-react with bovine, porcine rabbit and murine serum albumin.		
Concentration:	1 mg/ml		
UniProt ID:	P02768		
Specificity:	Human serum albumin		
Host:	Mouse		
Isotype:	Murine IgG1		
Clone:	AL-01		
Formulation:	Ascites diluted with PBS containing 0.1% sodium azide		
Applications:	WB, IHC, ELISA, IF		
Application Notes:	Suitable for affinity purification of albumin. Suitable for ELISA for detection of microalbuminuria		
Storage & Stability:	Store at 2-8°C for 1 month or in undiluted aliquots at -20°C to -80°C up to 1 year. Avoid repeated freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.