

AFP

Mouse Anti-Human Alpha-Fetoprotein (Clone AFP-01) mAb

| | | | |
|---------------------------------|---|------------------|--------|
| Catalog No. | MON5032A | Quantity: | 100 ug |
| Alternate Names: | Alpha-1-fetoprotein, Alpha-fetoglobulin | | |
| Description: | This antibody is highly specific for human alpha-fetoprotein. No reactivity with other serum proteins. The epitope is dependent on Ca ²⁺ and is sensitive to EDTA. | | |
| Concentration: | 0.5 mg/ml | | |
| UniProt ID: | P02771 | | |
| Specificity: | Human alpha-Fetoprotein | | |
| Host: | Mouse | | |
| Immunogen: | Purified alpha-Fetoprotein | | |
| Isotype: | IgG1 | | |
| Clone: | AFP-01 | | |
| Formulation: | PBS, pH 7.4 containing 15 mM sodium azide | | |
| Applications: | ELISA, IF, IP, WB | | |
| Storage & Stability: | Store at 2-8°C for 1 month or in undiluted aliquots at -20°C to -80°C up to 1 year. Avoid repeated freeze/thaw cycles. | | |

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

