

AFP

Mouse Anti-Human Alpha-Fetoprotein (Clone AFP-11) mAb

Catalog No.	MON5032B	Quantity:	100 ug
Alternate Names:	Alpha-1-fetoprotein, Alpha-fetoglobulin		
Description:	This antibody is highly specific for human alpha-fetoprotein. No reactivity with other serum proteins.		
Concentration:	0.5 mg/ml		
UniProt ID:	P02771		
Specificity:	Human alpha-Fetoprotein		
Host:	Mouse		
Immunogen:	Purified alpha-Fetoprotein		
Isotype:	IgG1		
Clone:	AFP-11		
Formulation:	PBS, pH 7.4 containing 15 mM sodium azide		
Applications:	ELISA, IF, IP, WB		
Application Notes:	Use for affinity purification. Use for detection in a sandwich ELISA with MON5032A.		
Storage & Stability:	Store at 2-8°C for 1 month or in undiluted aliquots at -20°C to -80°C up to 1 year. Avoid repeated freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

