

Mouse Anti-Horseradish Peroxidase (Clone HP-03) mAb

Catalog No.	MON5041	Quantity:	100 ug
Concentration:	1 mg/ml		
Specificity:	Horseradish peroxidase		
Immunogen:	Horseradish peroxidase, highly purified		
Isotype:	Murine IgG1		
Clone:	HP-03		
Formulation:	Purified antibody in PBS containing sodium azide		
Applications:	ELISA, ICC, IHC, WB		
Storage & Stability:	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.