

## Mouse Anti-Biotin (Clone BIO-08) mAb

<b>Catalog No.</b>	MON5042	<b>Quantity:</b>	1.0 mL
<b>Description:</b>	This antibody is useful for molecular biological applications and immunohistology. The antibody permits the formation of antibody-biotin-streptavidin-enzyme complexes thus enhancing the sensitivity of the detection system.		
<b>Concentration:</b>	0.1 mg/ml		
<b>Specificity:</b>	The antibody recognizes an epitope that is exposed by protein-biotin and nucleic acid-biotin conjugates.		
<b>Host:</b>	Mouse		
<b>Isotype:</b>	Murine IgG3		
<b>Clone:</b>	Bio-08		
<b>Formulation:</b>	Purified from tissue culture supernatant. Contains sodium azide.		
<b>Purification:</b>	Protein A affinity chromatography		
<b>Applications:</b>	IHC		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

