

## Mouse Anti-Streptavidin (Clone Strep-10) mAb

**Catalog No.** MON5043 **Quantity:** 1.0 mL

**Description:** The antibody can be used in molecular biological tests and immunohistology. The antibody permits the formation of antibody streptavidin-biotin-enzyme complexes thus enhancing the sensitivity of the detection system

**Concentration:** 0.25 mg/ml

**Specificity:** Streptavidin, does not react with avidin

**Isotype:** Murine IgG1

**Clone:** Strep-10

**Formulation:** Purified antibody solution containing sodium azide.

**Applications:** ELISA, IHC

**Storage & Stability:** Store at 2-8°C for up to 1 year.

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**