

## Rat Anti-Mouse IgA (Clone LO-MA-7) mAb

<b>Catalog No.</b>	MON5050	<b>Quantity:</b>	1 ml
<b>Description:</b>	This antibody is useful as second antibody in immunoassays, such as ELISA's and immunohistochemistry on frozen sections.		
<b>Concentration:</b>	1 mg/ml		
<b>Specificity:</b>	Mouse alpha heavy chain of immunoglobulin $K_a = 4.8 \times 10^9 \text{ mol}^{-1}$		
<b>Isotype:</b>	Rat IgM kappa		
<b>Clone:</b>	LO-MA-7		
<b>Formulation:</b>	Purified antibody solution in PBS with 0.1% sodium azide		
<b>Applications:</b>	IHC-F, ELISA		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.