

Rat Anti-Mouse IgA (Clone LO-MA-7) mAb, HRP

Catalog No.	MON5050P	Quantity:	1 ml
Description:	This antibody is useful as second antibody in immunoassays, such as ELISA's and immunohistochemistry on frozen sections.		
Concentration:	1 mg/ml		
Specificity:	Mouse alpha heavy chain of immunoglobulin $K_a = 4.8 \times 10^9 \text{ mol}^{-1}$		
Isotype:	Rat IgM kappa		
Clone:	LO-MA-7		
Conjugate:	Horseradish peroxidase (HRP)		
Formulation:	Purified antibody solution in PBS with 50% glycerol		
Applications:	IHC-F, ELISA		
Storage & Stability:	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

