

## Rat Anti-Mouse IgE (Clone LO-ME-3) mAb, Biotin

<b>Catalog No.</b>	MON5051B	<b>Quantity:</b>	1.0 mL
<b>Description:</b>	This antibody may be used as second antibody in immunoassays, such as ELISA and immunohistochemistry on frozen sections.		
<b>Concentration:</b>	1 mg/ml		
<b>Specificity:</b>	Mouse epsilon heavy chain of immunoglobulin $K_a = 3 \times 10^9 \text{ M}^{-1}$		
<b>Isotype:</b>	Rat IgG1 kappa		
<b>Clone:</b>	LO-ME-3		
<b>Conjugate:</b>	Biotin		
<b>Formulation:</b>	Purified antibody solution in PBS with 0.1% sodium azide		
<b>Applications:</b>	ELISA, IHC		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

