

Rat Anti-Mouse IgE (Clone LO-ME-3) mAb, HRP

Catalog No.	MON5051P	Quantity:	1.0 mL
Description:	This antibody may be used as second antibody in immunoassays, such as ELISA and immunohistochemistry on frozen sections.		
Concentration:	1 mg/ml		
Specificity:	Mouse epsilon heavy chain of immunoglobulin $K_a = 3 \times 10^9 \text{ M}^{-1}$		
Isotype:	Rat IgG1 kappa		
Clone:	LO-ME-3		
Conjugate:	Horseradish peroxidase		
Formulation:	Purified antibody solution in PBS with 50% glycerol		
Applications:	ELISA, IHC		
Storage & Stability:	Store at 2-8°C for up to 1 year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

