

Rat Anti-Mouse IgM (Clone LO-MM-9) mAb, FITC

Catalog No.	MON5057F	Quantity:	1.0 mL
Description:	This antibody may be used as second antibody in immunoassays, such as ELISA. Also used for IHC-F and immunoaffinity chromatography.		
Concentration:	1 mg/ml		
Specificity:	Mouse mu heavy chain of Immunoglobulin $K_a = 7 \times 10^8 \text{ M}^{-1}$		
Isotype:	Rat IgG2a kappa		
Clone:	LO-MM-9		
Conjugate:	FITC		
Formulation:	Purified antibody solution in PBS with 0.1% sodium azide		
Applications:	ELISA, IHC-F		
Storage & Stability:	Store at 2-8°C for up to 1 year. Protect from light.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

