

## Dinitrophenyl hapten, clone LO-DNP-2 Monoclonal Antibody

**Catalog No.:** MON 5070

**Quantity:** 1 ml

**Specificity**

Rat monoclonal antibody against Dinitrophenyl hapten.

Immunogen: DNP-ASCARIS, Nippostrongylus brasiliensis-DNP

Avidity:  $1.7 \times 10^{10} \text{ M}^{-1}$

**Immunoglobulin type**

Rat kappa IgG1, allotype IgK-1a

**Use**

ELISA

**Presentation**

1 ml with 1 mg purified antibody solution Catalog No.: MON 5050, (unconjugated).

FOR RESEARCH USE ONLY, NOT FOR DRUG, DIAGNOSTIC OR OTHER USE.



**Cell Sciences, Inc.**  
480 Neponset Street  
Bldg 12A  
Canton, MA 02021

Toll Free: 888-769-1246  
Phone: 781-828-0610  
Fax: 781-828-0542

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Web Site: [www.cellsciences.com](http://www.cellsciences.com)