

INTERLEUKIN-1 RECEPTOR TYPE 1 (p80), clone D14e3 Monoclonal Antibody

Catalog No.: MON 5078

Quantity: 1 ml

Specificity

This monoclonal antibody binds to soluble mouse interleukin-1 receptor type I. The IL-1 system includes two agonists (IL-1 α and IL-1 β), converting enzymes, antagonists, two receptors (IL-1 RI and IL-1 RII) and the IL-1 receptor accessory protein. Interleukin-1 signal is transduced through the type I receptor.

Immunoglobulin type

rat IgG1

Use

The anti mouse antibody can be used for Western Blots.

Presentation

1 ml (100 μ g/ml) purified antibody solution in PBS, 0.1% BSA, 0.02% Na-azide.

Literature

- Kojouharoff, G. et al., *Clin. Exp. Immunol.*, 1997, 107: 353-358.
- Stern, A. et al., *Arch. Biochem. Biophys*, 1989, 274(1): 26-36.

FOR RESEARCH USE ONLY, NOT FOR DRUG, DIAGNOSTIC OR OTHER USE.



Cell Sciences, Inc.
480 Neponset Street
Bldg 12A
Canton, MA 02021

Toll Free: 888-769-1246
Phone: 781-828-0610
Fax: 781-828-0542

E-mail: info@cellsciences.com
Web Site: www.cellsciences.com