

Rat Anti-Mouse Ig Gamma light chain (Clone LO-MK-1) mAb

Catalog No.	MON5083
Lot No:	1209901
Quantity:	1 ml
Specificity:	Mouse gamma light chain of Immunoglobulin (determined by immunodot) and exhibits no reactivity to mouse lambda light chain. No detectable reactivity with human IgG as tested by ELISA.
Concentration:	1.0 mg/ml purified antibody
Host:	Rat
Isotype:	Rat IgG2a, kappa
Clone:	LO-MK-1
Conjugate:	Unconjugated MON5083B with Biotin label MON5083P with HRP label
Formulation:	PBS with 50% Glycerol
Applications:	Useful as second antibody in immunoassays, such as ELISA (good binding on plastic), immunohistochemistry, flow cytometry, immunofluorescence and Western blotting.
Storage & Stability:	Store antibody at 2-8°C until expiration date.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

