

Mouse Anti-Hepatitis B Virus Antigen Clone HB5 mAb

Catalog No.	MON6081	Quantity:	100 µg
Description:	Recognizes following Hepatitis B virus subtypes: ayw2, ayw3, ayw4, ayr, adw2, adw4, adrq+, adrq-. HB5 does not cross-block with HB3.		
Host:	Mouse		
Isotype:	IgG2a		
Clone:	HB5		
Formulation:	100 µg in solution (PBS containing 15 mM sodium azide, pH 7.4). Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
Applications:	The antibody can be used for ELISA.		
Storage & Stability:	Store at 4°C for 3 months. For long-term storage aliquot and store at -20°C. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

