

## EP84R CDS

### Mouse Anti-ASFV EP84R (1A025) mAb

<b>Catalog No.</b>	MVV11202A MVV11202B	<b>Quantity:</b>	50 µg 100 µg
<b>Alternate Names:</b>	EP84R, EP84R protein, Protein EP84R, African swine fever virus (ASFV)		
<b>Description:</b>	Mouse Anti-ASFV EP84R (1A025) is a murine monoclonal antibody which targets protein EP84R from African swine fever virus (ASFV)		
<b>Accession (target):</b>	YP_009927178.1		
<b>Source:</b>	Mammalian cells		
<b>Specificity:</b>	African swine fever virus (ASFV) protein hypothetical protein IM014_gp088, EP84R CDS		
<b>Immunogen:</b>	Recombinant ASFV EP84R protein		
<b>Brand:</b>	Antibody System		
<b>Isotype:</b>	Mouse IgG2a		
<b>Formulation:</b>	Sterile filtered 0.01M PBS buffer, pH 7.4 containing 0.05% Proclin 300, 50% glycerol		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Purification:</b>	Protein A/G affinity chromatography from cell culture supernatant		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Clone:</b>	1A025		
<b>Applications:</b>	ELISA, Western blot		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences®**  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)