

Casp12

Synthetic Mouse Caspase 12 (aa 2-17)(NT) Blocking Peptide

| | | | |
|---------------------------------|---|------------------|-------|
| Catalog No. | PX037BP | Quantity: | 50 µg |
| Description: | Amino acids 2 to 17 of mouse caspase-12. | | |
| Gene ID: | 12364 | | |
| Application: | The peptide is used for blocking the antibody activity of caspase-12. Incubating the peptide with equal volume of antibody for 30 min at 37°C usually completely blocks the antibody activity in Western blotting. | | |
| Formulation: | It is supplied as 200 µg/ml, 50 µg/vial , in PBS pH7.2 (10 mM NaH ₂ PO ₄ , 10 mM, Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide.. Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only. | | |
| Sequence: | NH ₂ -AARRTHERDPIYKIKG-OH | | |
| Storage & Stability: | Store at -20°C, stable for one year. | | |

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

