


## Ripk3

### Synthetic Mouse Receptor-Interacting Protein Kinase 3 (aa 473-486) Blocking Peptide

<b>Catalog No.</b>	PX096BP	<b>Quantity:</b>	50 µg
<b>Alternate Names:</b>	2610528K09Rik, AW107945, Rip3, receptor interacting protein 3		
<b>Description:</b>	Amino acids 473 to 486 of murine RIP3.		
<b>Gene ID:</b>	56532		
<b>Application:</b>	The peptide is used for blocking the antibody activity of RIP3. The peptide with equal volume of antibody for 30 min at 37°C usually completely blocks the antibody activity in Western blotting.		
<b>Formulation:</b>	It is supplied as 200 µg/ml, 50 µg/vial , in PBS pH7.2 (10 mM NaH <sub>2</sub> PO <sub>4</sub> , 10 mM, Na <sub>2</sub> HPO <sub>4</sub> , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide.. <b>Precaution:</b> Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
<b>Sequence:</b>	NH <sub>2</sub> -AQFGRGRGWQPFHK-OH		
<b>Storage &amp; Stability:</b>	Store at -20°C, stable for one year.		

RIP3 Peptide



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

