

RYBP

Synthetic RYBP/DEDAF (aa 215-228)

Catalog No.	PX202BP	Quantity:	50 µg
Alternate Names:	AAP1, DEDAF, YEAF1, YY1 and E4TF1 associated factor 1, apoptin-associating protein 1, death effector domain-associated factor, ring1 interactor RYBP		
Description:	Amino acids 215 to 228 of human DEDAF.		
Gene ID:	23429		
Application:	The peptide is used for blocking the activity of DEDAF. The peptide with equal volume of antibody for 30 min at 37°C usually completely blocks the antibody activity in Western blotting.		
Formulation:	It is supplied as 200 µg/ml, 50 µg/vial , in PBS pH7.2 (10 mM NaH ₂ PO ₄ , 10 mM, Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide.. Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
Sequence:	NH ₂ -TPKGDMSAVNDESF-OH		
Storage & Stability:	Store at -20°C, stable for one year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

