

Streptavidin - HRP Conjugate

Catalog No.	PX741	Quantity:	0.5 mg
Alternate Names:	Streptavidin-horseradish peroxidase		
Description:	Streptavidin-HRP conjugate to be used for ELISA and Western Blot analysis		
Physical State:	Liquid		
Formulation:	0.01 M Na ₂ PO ₄ , 0.25 M NaCl, pH 7.6 with 25% glycerol		
Applications:	Use at 1:2,000 - 1:10,000 dilution for ELISA or Western blot. Remaining working strength solution should be discarded at the end of the day. Optimum dilution should be determined for each application by the end-user.		
Storage & Stability:	Upon receipt, store at 2-8° C for up to 2 weeks. It is recommended to store as single-use aliquots (≥ 20 µl) at -80° C.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.