

SERPINB14

Human Anti-OVA/Ovalbumin (SAA0619) mAb

Catalog No.	RCK15101A RCK15101B	Quantity:	100 µg 1.0 mg
Alternate Names:	Ovalbumin, Allergen Gal d II, Egg albumin, Plakalbumin		
Description:	Human Anti-OVA/Ovalbumin (SAA0619) mAb is purified from cell culture supernatant by affinity chromatography.		
UniProt ID (target):	P01012		
Target:	Ovalbumin from Chicken <i>Gallus gallus</i>		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Human IgG1		
Formulation:	Sterile filtered 0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Clone:	SAA0619		
Purification:	Protein A/G affinity chromatography from cell culture supernatant.		
Applications:	ELISA, Western blot EC ₅₀ = 5.956 ng/ml, clone SAA0619 binds with ovalbumin by indirect ELISA		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com