

Anti-Fel d 1 Human IgE (1B7) mAb

Catalog No.	RFK29907A RFK29907B	Quantity:	100 µg 1.0 mg
Alternate Names:	Fel d 1, Felis domesticus allergen		
Description:	Affinity purified human monoclonal IgE antibody (clone 1B7) that targets Fel d 1, Felis domesticus allergen from Felis catus (Cat) (Felis silvestris catus)		
Target:	Fel d 1, Felis domesticus allergen from Felis catus (Cat) (Felis silvestris catus)		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Human IgE		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Clone:	1B7		
Applications:	ELISA		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

