

Mouse Anti-BC2 tag (PDRKAAVSHWQQ) (SAA2177) mAb

Catalog No.	RGE07602A RGE07602B	Quantity:	50 µg 100 µg
Alternate Names:	BC2-tag, PDRKAAVSHWQQ, BC2 tag		
Description:	Mouse monoclonal antibody which targets the BC2 tag sequence PDRKAAVSHWQQ		
UniProt ID:	P35222		
Specificity:	BC2 tag sequence PDRKAAVSHWQQ		
Concentration:	1.0 mg/ml, lot specific		
Isotype:	Mouse VHH-mFc		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Source:	Mammalian cell culture		
Purity:	> 95% by SDS-PAGE.		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	ELISA		
Storage & Stability:	Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com