

Rabbit Anti-6X His Tag (3D5) mAb

Catalog No.	RGK07901B RGK07901C	Quantity:	100 µg 1.0 mg
Alternate Names:	6X His Tag, 6His 6xHis, HHHHHH		
Description:	Rabbit monoclonal anti-6X His Tag (3D5) was expressed in mammalian cell line and purified using Protein A or G affinity chromatography.		
PubChem:	65396		
Concentration:	1.0 mg/ml, lot specific		
Isotype:	Rabbit IgG1 kappa		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Clone:	3D5		
Purity:	> 95% by SDS-PAGE.		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	Western blot, IF, IP		
Storage & Stability:	Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

