

## Mouse Anti-PEG (2B5) mAb

<b>Catalog No.</b>	RGK24706A	<b>Quantity:</b>	100 µg
	RGK24706B		1.0 mg

**Alternate Names:** PEG, Poly(ethylene glycol), Polyethylene glycol

**Description:** Polyethylene glycol (PEG) is a polymer which is widely used for covalent modification of biological polymers such as proteins and peptides. PEG is often used in drug delivery and drug modification technologies, which can be directly coupled to the drug, or attached to the drug surface and encapsulated in nanomaterials. A full understanding of the pharmacokinetic (PK) characteristics and biological analysis of PEG and PEG drugs is of great significance to evaluate their overall safety and promote the development of related technology. PEG antibody is an important tool to detect PEG-modified drugs, in PK studies and immunogenicity analysis.

**Concentration:** 1.0 mg/ml

**Isotype:** Mouse IgG1 kappa

**Formulation:** Sterile-filtered 0.01M PBS, pH 7.4

**Clone:** 2B5

**Purity:** > 95% by SDS-PAGE.

**Purification:** Protein A/G affinity chromatography from cell culture supernatant

**Applications:** ELISA, FC

**Storage & Stability:** Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. **Avoid repeated freeze-thaw cycles.**

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences**<sup>®</sup>  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)