

Human (IgE) Anti-Alpha-gal (16D9) mAb

Catalog No.	RGK25408A RGK25408B	Quantity:	100 µg 1.0 mg
Alternate Names:	Galactose-α-1,3-galactose		
Description:	Human monoclonal antibody, isotype IgE, that targets Alpha-gal, was purified by affinity chromatography.		
Specificity:	Alpha-gal, general specificity		
Concentration:	1.0 mg/ml, lot specific		
Isotype:	Human IgE		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Clone:	16D9		
Purity:	> 95% by SDS-PAGE.		
Purification:	Affinity chromatography from cell culture supernatant		
Applications:	ELISA		
Storage & Stability:	Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

