

Mouse Anti-polyglutaminepolyQ (MW1) mAb

Catalog No.	RGK26001A RGK26001B	Quantity:	100 µg 1.0 mg
Alternate Names:	Glutamine, Poly-Glutamine, Polyglutamine, poly-Q, polyQ, GQQQQQQQQQQG		
Description:	Highly specific, affinity purified, mouse monoclonal antibody that targets expanded PolyQ repeats from Huntingtin (Htt)		
ChEBI:	28300		
Specificity:	Polyglutamine, GQQQQQQQQQQG		
Concentration:	1.0 mg/ml, lot specific		
Isotype:	Mouse IgG1		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Clone:	MW1		
Purity:	> 95% by SDS-PAGE.		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	ELISA, ICC, IP, Western blot		
Storage & Stability:	Store at 2-8°C for 1-2 weeks or unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com