

Mouse Anti-Twin-Strep-tag (SAA0348) mAb

Catalog No.	RGK26103A RGK26103B	Quantity:	50 µg 100 µg
Alternate Names:	Twin Strep tag, SAWSHQPFEKGGGSGGGSGGSAWSHQPFEK		
Description:	Mouse monoclonal antibody which targets the Twin-Strep tag.		
Specificity:	Twin-Strep tag		
Concentration:	1.0 mg/ml, lot specific		
Isotype:	Mouse IgG2b		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Source:	Mammalian cell culture		
Purity:	> 95% by SDS-PAGE.		
Purification:	Protein A/G affinity chromatography		
Applications:	ELISA, WB		
Storage & Stability:	Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.