

## Mouse Anti-ALFA Tag (SAA2164) mAb

<b>Catalog No.</b>	RGK26502A RGK26502B	<b>Quantity:</b>	50 µg 100 µg
<b>Alternate Names:</b>	ALFA-tag, (P)SRLEEEELRRRLTE(P), SRLEEEELRRRLTE, ALFA tag		
<b>Description:</b>	VHH-mFc monoclonal antibody which targets the ALFA tag sequence.		
<b>Specificity:</b>	ALFA tag sequence		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Isotype:</b>	VHH-mFc		
<b>Formulation:</b>	Sterile-filtered 0.01M PBS, pH 7.4		
<b>Source:</b>	Mammalian cell culture		
<b>Purity:</b>	> 95% by SDS-PAGE.		
<b>Purification:</b>	Protein A/G affinity chromatography		
<b>Clone:</b>	SAA2164		
<b>Applications:</b>	ELISA, FCM, IHC, IP, WB		
<b>Storage &amp; Stability:</b>	Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

