

Mouse Anti-Pan-pS(pSer)/pT(pThr) (SAA0308) mAb

Catalog No.	RGK26902A RGK26902B	Quantity:	100 µg 1.0 mg
Alternate Names:	Phospho-specific peptide, Pan-pS(pSer)/pT(pThr)		
Description:	Mouse monoclonal antibody which targets a broad spectrum of phospho-specific peptides, Pan-pS(pSer)/pT(pThr)		
Specificity:	Phospho-specific peptides, pS(pSer)/pT(pThr)		
Concentration:	1.0 mg/ml, lot specific		
Isotype:	Mouse IgG2a, kappa		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Source:	Mammalian cell culture		
Purity:	> 95% by SDS-PAGE.		
Purification:	Protein A/G affinity chromatography		
Clone:	SAA0308		
Applications:	ELISA, WB, SPR		
Storage & Stability:	Store unopened at 2-8°C for 1 week or -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

