



Mouse Anti-EPEA Tag / C-tag (SAA0356) mAb

Catalog No.	RGK30901A RGK30901B	Quantity:	50 µg 100 µg
Alternate Names:	EPEA Tag (GGEPEA), EPEA Tag, C-tag		
Description:	Mouse monoclonal antibody which targets the EPEA-tag (also known as C-Tag) which can only be used at the C-terminus of target proteins.		
Specificity:	EPEA Tag (GGEPEA)		
Concentration:	1.0 mg/ml, lot specific		
Isotype:	VHH-mFc		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Source:	Mammalian cell culture		
Purity:	> 95% by SDS-PAGE and SEC-HPLC		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Clone:	SAA0356		
Applications:	ELISA, WB		
Manufacturer:	Antibody System		
Storage & Stability:	Store unopened at 2-8°C for 1 week or -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		



NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.