

Mouse Anti-STX/Saxitoxin (5F7) mAb

Catalog No.	RGK32902A RGK32902B	Quantity:	100 µg 1.0 mg
Alternate Names:	Saxitoxin, STX		
Description:	Mouse anti-Saxitoxin monoclonal antibody is purified by Protein A or G affinity chromatography from cell culture supernatant.		
Accession:	ChEBI: 34970		
Concentration:	1.0 mg/ml, lot specific		
Isotype:	Mouse IgG2a kappa		
Formulation:	Sterile-filtered 0.01M PBS, pH 7.4		
Clone:	5F7		
Purity:	> 95% by SDS-PAGE.		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	ELISA, SPR		
Storage & Stability:	Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

