

Mouse Anti-G-quadruplex DNA (8H2) mAb

Catalog No.	RGK60302A	Quantity:	100 µg
	RGK60302B		1.0 mg

Alternate Names: DNA, G-quadruplex DNA

Description: Anti-DNA/RNA G-quadruplex structures, clone BG4 antibody, is a highly specific single-chain antibody that targets DNA/RNA G-quadruplex structures and has been tested in EMSA, ELISA, IF and ChIP-seq

This recombinant antibody has high specificity, high affinity, and high reproducibility without hybridoma cells, preservatives or animal components.

Concentration: 1.0 mg/ml

Isotype: Mouse IgG2a kappa

Formulation: Sterile-filtered 0.01M PBS, pH 7.4

Clone: 8H2

Purity: > 95% by SDS-PAGE.

Purification: Protein A/G affinity chromatography from cell culture supernatant

Applications: ELISA, IF, Immunoblot

Storage & Stability: Store unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. **Avoid repeated freeze-thaw cycles.**

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com