

## Mouse Anti-DA/Domoic acid (24cB7) mAb

<b>Catalog No.</b>	RGK81801A RGK81801B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Domoic acid, DA, DOQ, L-domoic acid, CAS:14277-97-5		
<b>Description:</b>	Monoclonal mouse anti-DA/Domoic acid (24cB7) was affinity purified from cell culture supernatant.		
<b>ChEBI:</b>	34727		
<b>Specificity:</b>	Domoic acid		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Isotype:</b>	Mouse IgG2a kappa		
<b>Formulation:</b>	Sterile-filtered 0.01M PBS, pH 7.4		
<b>Clone:</b>	24cB7		
<b>Purity:</b>	> 95% by SDS-PAGE and SEC-HPLC		
<b>Purification:</b>	Protein A/G affinity chromatography from cell culture supernatant		
<b>Applications:</b>	ELISA		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for 1-2 weeks or unopened at -20°C to -80°C for up to 1 year. It is recommended to prepare aliquots and store at -20°C to -80°C. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences**<sup>®</sup>  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)